Homedics

Heated Shiatsu Back Massager

Register your product today at
www.homedics.co.uk/product-registration

3 YEAR GUARANTEE

SBM-150H-WW

**READ ALL INSTRUCTIONS BEFORE USE.
SAVE THESE INSTRUCTIONS FOR FUTURE
REFERENCE.**

- This appliance is suitable for ages 8 and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- ALWAYS unplug the USB cable from the electrical outlet immediately before cleaning. To turn off, press and hold ⏻.
- NEVER leave an appliance unattended when plugged in. Unplug from the mains outlet when not in use and before putting on or taking off parts or attachments.
- Do not place or store appliance where it can fall or be pulled into a bath or sink. Do not place in or drop into water or other liquid.
- DO NOT reach for an appliance that has fallen in water or other liquids. Switch off at the mains and unplug it immediately. Keep dry – DO NOT operate in wet or moist conditions.
- NEVER insert pins, metallic fasteners or objects into the appliance or any opening.
- Use this appliance for the intended use as described in this booklet. DO NOT use attachments not recommended by Homedics.
- For indoor use only.
- This product requires a 5V DC 2A USB Certified power supply.
- NEVER operate the appliance if charging cable is damaged, if it is not working properly, if it has been dropped or damaged, dropped into water.
- Keep cord away from heated surfaces.
- DO NOT operate where aerosol (spray) products are being used or where oxygen is being administered.
- DO NOT operate under a blanket or pillow. Excessive heating can occur and cause fire, electrocution or injury to persons.
- DO NOT carry this appliance by the cord or use cord as a handle.
- Use the carry handle supplied on the product to carry.
- DO NOT attempt to repair the appliance. There are no user serviceable parts. All servicing of this appliance must be performed at an authorised Homedics Service Centre.
- NEVER sit or stand on the massage part (back) of the device. The unit should only be used when attached to a chair with the back in the vertical position.
- Never block the air openings of the appliance. Keep air openings free of lint, hair etc.
- Burns can result from improper use.
- If applicable, the back strap should only be used to hold the cushion against the back of a chair. It is not intended for any other use.
- Please ensure that all hair, clothing and jewellery are kept clear of the massage mechanism or any other moving parts of the product at all times.
- If you have any concerns regarding your health, consult a doctor before using this appliance.
- This is a non-professional appliance, designed for personal use and intended to provide a soothing massage to worn muscles. DO NOT use as a substitute for medical attention.

- Use of this product should be pleasant and comfortable. Should pain or discomfort result, discontinue use and consult your GP.
- Pregnant women and diabetics should consult a doctor before using this item. Not recommended for use by individuals with sensory deficiencies including diabetic neuropathy.
- DO NOT use on an infant, invalid or on a sleeping or unconscious person.
- DO NOT use on insensitive skin or on a person with poor blood circulation.
- This appliance should NEVER be used by any individual suffering from any physical ailment that would limit the user's capacity to operate the controls.
- NEVER use directly on swollen or inflamed areas or skin eruptions.
- Do not use on irritated, sunburned or chapped skin, or on any area with dermatological problems.
- DO NOT use before going to bed. The massage has a stimulating effect and can delay sleep.
- Do not use for longer than the recommended time.
- Only gentle force should be exerted against the massage mechanism in order to eliminate risk of injury.
- The appliance has a heated surface. Persons insensitive to heat must be careful when using the appliance.
- If the product feels excessively hot while charging, switch off and unplug the USB cable from the electrical outlet immediately and contact the Homedics Service Centre.
- Excessive use could lead to the product's excessive heating and shorter life. Should this occur, discontinue use and allow the unit to cool before operating.
- This product cannot be used while charging.
- Do not crease the appliance by placing items on top of it during storage.
- This appliance contains batteries that are non-replaceable.
- The battery must be removed from the appliance before it is scrapped.
- The battery is to be disposed of safely.
- Batteries (battery pack or batteries installed) shall not be exposed to excessive heat such as sunshine, fire or the like.
- Failure to follow the above may result in the risk of fire or injury.

3 YEAR GUARANTEE

FKA Brands Ltd guarantees this product from defect in material and workmanship for a period of 3 years from the date of purchase, except as noted below. This FKA Brands Ltd product guarantee does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorised accessory; alteration to the product; or any other conditions whatsoever that are beyond the control of FKA Brands Ltd. This guarantee is effective only if the product is purchased and operated in the UK / EU. A product that requires modification or adaptation to enable it to operate in any country other than the country for which it was designed, manufactured, approved and / or authorised, or repair of products damaged by these modifications is not covered under this guarantee. FKA Brands Ltd shall not be responsible for any type of incidental, consequential or special damages.

To obtain guarantee service on your product, return the product post-paid to your local service centre along with your dated sales receipt (as proof of purchase). Upon receipt, FKA Brands Ltd will repair or replace, as appropriate, your product and return it to you, post-paid. Guarantee is solely through Homedics Service Centre. Service of this product by anyone other than the Homedics Service Centre voids the guarantee. This guarantee does not affect your statutory rights. For your local Homedics Service Centre, go to www.Homedics.co.uk/servicecentres

WEEE explanation

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

Battery directive

This symbol indicates that batteries must not be disposed of in the domestic waste as they contain substances which can be damaging to the environment and health. Please dispose of batteries in designated collection points.

Battery replacement

Your product includes a rechargeable battery designed to last the lifetime of the product. In the unlikely event that you should require a replacement battery, please contact Customer Services, who will supply details of the warranty and out-of-warranty replacement battery service.

PRODUCT FEATURES

- 1. Carry Handle
- 2. Massage cushion
- 3. Elastic strap
- 4. Lower back warming pad
- 5. Controller
- 6. (M) Mode button
- 7. (W) Warming button
- 8. (P) Power button
- 9. USB charging LED
- 10. USB charging port

CHARGING

- Please fully charge the product before first use.
- Connect the USB cable to the charging port and the other end to a USB power source (5V/2A).
- The charging LED will flash while charging and remain solid when fully charged. A full charge will take approximately 5 hours.
- Run time from a full charge is up to 3 hrs 20 mins without heat, and up to 2hrs 20 mins with heat.

NOTE: The device cannot be used while charging.

INSTRUCTIONS FOR USE

- Place the massager on a chair or sofa, using the elastic strap to hold it in place if required.
- Sit with your back against the massager, ensuring you are centrally aligned.
- Press the power button to start the massage experience. (Power LED = white)
- Press the power button again to increase the massage speed. (Power LED = red)
- Press the (M) button to scroll through the four massage programs.

1	mode LEDs both lit	Full back massage
2	mode LEDs both flashing	Spot massage
3	upper mode LED lit	Upper back massage
4	lower mode LED lit	Lower back massage

- The warming function for the lower back area is automatically activated when the product is first turned on and the red warmth LED will be lit. To turn the warming feature off press the warming button (LED turns off). Press again to turn back on.
- Each program runs for 15 minutes. Press and hold the power button to turn off sooner.

HINTS AND TIPS

- If you prefer a gentler massage, placing a folded towel between your back and the massager will soften the massage pressure.
- Moving the massager higher up, or lower down on your chair allows you to easily target different areas of the back
- The warming pad in the lower back area is designed to be used alongside your massage. It is not possible to use the warming function without massage.

CLEANING AND CARE

- Unplug the appliance and allow to cool before cleaning.
- Clean only with a soft, slightly damp cloth.
- Never allow water or any other liquids to come into contact with the unit.
- DO NOT immerse in any liquid to clean.
- NEVER use abrasive cleaners, brushes, or strong chemicals.
- Store in a cool dry place when not in use.

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FKA Brands Ltd, Somerhill Business Park, Tonbridge, Kent TN11 0GP, UK
FKA Brands Ltd, 29 Earlsfort Terrace, Dublin D02AY28, Ireland